



LiNA PORT

Single use 5mm
Cannula and Trocar System

Ordering Information:

CE
LiNA Port 0459

REF.	Product: Dilating Tip	Packing
1005-06D1	5mm Cannula and Dilating Trocar - short	Single pack, sterile / Box of 5 sets
1005-06D2	5mm Cannula and Dilating Trocar - short	Dual pack*, sterile / Box of 5 sets
1005-09D1	5mm Cannula and Dilating Trocar - long	Single pack, sterile / Box of 5 sets
1005-09D2	5mm Cannula and Dilating Trocar - long	Dual pack*, sterile / Box of 5 sets
1005-06D1-GV	5mm Cannula with gas port and on/off valve and Dilating Trocar - short	Single pack, sterile / Box of 5 sets
1005-06D2-GV	5mm Cannula with gas port and on/off valve + one Cannula without gas port + Dilating Trocar - short	Dual pack*, sterile / Box of 5 sets
1005-09D1-GV	5mm Cannula with gas port and on/off valve and Dilating Trocar - long	Single pack, sterile / Box of 5 sets
1005-09D2-GV	5mm Cannula with gas port and on/off valve + one Cannula without gas port + Dilating Trocar - long	Dual pack*, sterile / Box of 5 sets
1005-06D1-G	5mm Cannula with gas port without on/off valve and Dilating Trocar - short	Single pack, sterile / Box of 5 sets
1005-06D2-G	5mm Cannula with gas port without on/off valve + one Cannula without gas port + Dilating Trocar - short	Dual pack*, sterile / Box of 5 sets
1005-09D1-G	5mm Cannula with gas port without on/off valve and Dilating Trocar - long	Single pack, sterile / Box of 5 sets
1005-09D2-G	5mm Cannula with gas port without on/off valve + one Cannula without gas port + Dilating Trocar - long	Dual pack*, sterile / Box of 5 sets

REF.	Product: Pyramidal Tip	Packing
1005-1	5mm Cannula and Pyramidal Trocar - short	Single pack, sterile / Box of 5 sets
1005-2	5mm Cannula and Pyramidal Trocar - short	Dual pack*, sterile / Box of 5 sets
2005-1	5mm Cannula and Pyramidal Trocar - long	Single pack, sterile / Box of 5 sets
2005-2	5mm Cannula and Pyramidal Trocar - long	Dual pack*, sterile / Box of 5 sets
1008-1	5mm Cannula with gas port and on/off valve and Pyramidal Trocar - short	Single pack, sterile / Box of 5 sets
1008-2	5mm Cannula with gas port and on/off valve + one Cannula without gas port + Pyramidal Trocar - short	Dual pack*, sterile / Box of 5 sets
2008-1	5mm Cannula with gas port and on/off valve and Pyramidal Trocar - long	Single pack, sterile / Box of 5 sets
2008-2	5mm Cannula with gas port and on/off valve + one Cannula without gas port + Pyramidal Trocar - long	Dual pack*, sterile / Box of 5 sets
1009-1	5mm Cannula with gas port without on/off valve and Pyramidal Trocar - short	Single pack, sterile / Box of 5 sets
1009-2	5mm Cannula with gas port without on/off valve + one Cannula without gas port + Pyramidal Trocar - short	Dual pack*, sterile / Box of 5 sets
2009-1	5mm Cannula with gas port without on/off valve and Pyramidal Trocar - long	Single pack, sterile / Box of 5 sets
2009-2	5mm Cannula with gas port without on/off valve + one Cannula without gas port + Pyramidal Trocar - long	Dual pack*, sterile / Box of 5 sets

REF.	Product: Conical Tip	Packing
1007-1	5mm Cannula and Conical Trocar - short	Single pack, sterile / Box of 5 sets
1007-2	5mm Cannula and Conical Trocar - short	Dual pack*, sterile / Box of 5 sets
2007-1	5mm Cannula and Conical Trocar - long	Single pack, sterile / Box of 5 sets
2007-2	5mm Cannula and Conical Trocar - long	Dual pack*, sterile / Box of 5 sets

*2 x Cannula + 1 x Trocar + 1 Clear Sight Wipe



- Lightweight cannula with or without gas port
- Choice of Pyramidal, Conical or Dilating Trocar tip
- Silicone valve - airtight seal
- Threaded cannula - efficient abdominal retention
- Reduced procedure cost - cost effective
- Sterile, single use - no cleaning, no risk of cross contamination

LiNA Single Use Cannula and Trocar

The LiNA Single Use Cannula and Trocar System meets the demand for quality and efficiency at the most competitive price. The small size allows easy access to the abdomen and easy manoeuvrability of instruments.

Two Lengths - Short and Long Ports: Short 62mm and Long 88mm. Diameter inside cannula: ø5,8mm allows easy insertion of all 5mm instruments. Suitable for all laparoscopic procedures and all patients.

Gas Ports: The cannulas are available with either a simple gas port, a gas port with on/off valve or without gas port.

Three Tips: The conical and pyramidal tip for sharp and easy insertion and the blunt dilating tip for safer entry.

Single Pack and Dual Pack: The LiNA Port is available in Single pack and Dual pack. The Dual Pack includes one trocar, two cannulas and one Clear Sight Wipe Anti-Fog Solution for camera optics.

Single Pack



Dual Pack with one LiNA Clear Sight Wipe

Airtight Silicone Double Seal Valve.

Secures position of instrument. Easy exchange of instruments. Enables the use of instruments down to 3mm.



Non-traumatic Fascia Thread.
Safe and firm anchoring of cannula.

Dilating, Conical or Pyramidal Trocar tip.
Supports the most used entry techniques.