



LiNA Surgipen

Reusable finger-controlled electrosurgical pencil for open monopolar electrosurgery

✓ **safe**

✓ **versatile**

✓ **functional**





functionality

LiNA Surgipen

The slim lightweight design makes it easy and comfortable for the surgeon to use during long procedures. A soft and flexible silicone cable ensures easy handling and long durability. The silicone cable remains soft even after many autoclavings.



safety

Press Button switches with "click"

Sealed waterproof switches ensure durability and safety in the operating theatre. The switches' homogenous activation pressure provides very high comfort and reduced movement of the pencil when activated.



versatility

Versatile Product Range

A wide range of possible combinations in connectors, adaptors, cable lengths and electrodes makes it possible to meet all operating theatre demands. The LiNA Surgipen can be used for all major electrosurgical units either with the 3-pin connector, ICC connector or by using adaptors.

Ordering information:

LiNA Surgipen **CE** 0459

REF:	Product	Packing
508-L-3	LiNA Surgipen - Reusable - Autoclavable - for 2.35mm electrodes. With 3m silicone cable and standard US-3-pin connector. With 70mm blade electrode.	Box of 10
508-L-5	LiNA Surgipen - Reusable - Autoclavable - for 2.35mm electrodes. With 5m silicone cable and standard US-3-pin connector. With 70mm blade electrode.	
508-L-3-ACC	LiNA Surgipen - Reusable - Autoclavable - for 2.35mm electrodes. With 3m silicone cable and ICC connector for Erbe ICC generators. With 70mm blade electrode.	Box of 10
508-L-5-ACC	LiNA Surgipen - Reusable - Autoclavable - for 2.35mm electrodes. With 5m silicone cable and ICC connector for Erbe ICC generators. With 70mm blade electrode.	
512-L-3	LiNA Surgipen - Reusable - Autoclavable - for 4mm hexagonal electrodes. With 3m silicone cable and standard US-3-pin connector.	Box of 10
512-L-5	LiNA Surgipen - Reusable - Autoclavable - for 4mm hexagonal electrodes. With 5m silicone cable and standard US-3-pin connector.	
512-L-3-ACC	LiNA Surgipen - Reusable - Autoclavable - for 4mm hexagonal electrodes. With 3m silicone cable and ICC connector for Erbe ICC generators.	Box of 10
512-L-5-ACC	LiNA Surgipen - Reusable - Autoclavable - for 4mm hexagonal electrodes. With 5m silicone cable and ICC connector for Erbe ICC generators.	
SET-0508	LiNA Smoke evacuation tubing set - Single use - 3m tubing with a blood/fluid stop filter and 22mmØ connector for LiNA, Erbe, Conmed, Stackhouse, Buffalo, Atmos, Finesse and Aaron smoke evacuation.	Box of 10 sterile
580	LiNA Tip Cleaner - Single use - Self Adhesive. 45x50mm with X-ray detected thread.	Box of 30 sterile
519	LiNA Pencil Holster - Semi reusable. Autoclavable.	Box of 10
571-ACC	LiNA HF generator Adaptor - reusable - For Erbe ICC generators.	Box of 1
571-AE	LiNA HF generator Adaptor - reusable - For Aesculab generators.	
571-B	LiNA HF generator Adaptor - reusable - For Bovie generators.	
571-E	LiNA HF generator Adaptor - reusable - For Old Erbe generators.	
571-L	LiNA HF generator Adaptor - reusable - For Lamidey generators.	
571-M	LiNA HF generator Adaptor - reusable - For Martin generators.	
571-N	LiNA HF generator Adaptor - reusable - For Neomed generators.	
571-S1	LiNA HF generator Adaptor - reusable - For Siemens type 1 generators.	
571-S2	LiNA HF generator Adaptor - reusable - For Siemens type 2 generators.	
572	LiNA HF generator Adaptor - reusable 8mm connector for foot controlled pencils - REF 619/517.	

